
ACADEMIE IT EN MEDIADESIGN

# EYE TRACKING MET TOBII PRO LAB

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# TEST OPNEMEN - CALIBRATIE

The screenshot shows a software interface for calibration. On the left, there is a dark area with several white crosshair markers. On the right, there is a panel with the following sections:


- View options:** Two radio buttons, 'Average results' (selected) and 'Left and right eyes results'.
- Calibration average results:** A table with columns for Calibration Date, Accuracy, Precision (SD), Precision (RMS), and Data loss.
- Buttons:** 'Recalibrate all points', 'Revalidate', and 'Accept'.

Orange callout boxes highlight the buttons: 'Recalibrate all points' (pointing to the button), 'Accept the results' (pointing to the 'Accept' button), and 'Revalidate specific points' (pointing to the 'Revalidate' button).

Calibration Date	Accuracy			Precision (SD)			Precision (RMS)			Data loss percent
	deg	px	mm	deg	px	mm	deg	px	mm	
11/3/2020 4:00 PM	0.05°	2	0.6	0.10°	4	1.2	0.14°	6	1.6	0%
11/3/2020 1:41 PM	0.06°	3	0.7	0.14°	6	1.7	0.17°	7	2.0	0%
11/3/2020 1:41 PM	0.10°	4	1.2	0.24°	10	2.8	0.24°	10	2.8	0%
11/3/2020 1:40 PM	0.05°	2	0.6	0.09°	4	1.1	0.11°	5	1.4	0%
11/3/2020 1:38 PM	0.14°	3	1.5	0.20°	4	1.1	0.20°	4	1.1	0%

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
# OPNAME STOPPEN

**i** Recording finished

Recording1 with Participant Participant1 was finished successfully.

You have to get participant consent to store and/or transfer the recorded eye tracking data. More information at:

<https://transparency.tobii.com/> 

- Save recording**  
The participant gave consent to process their eye tracking data. The recording and participant information will be saved.
- Delete recording**  
The participant did not give consent to process their eye tracking data. The recording and participant information will be deleted.

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